


# Self-supervised Denoising and Bulk Motion Artifact Removal of 3D Optical Coherence Tomography Angiography of Awake Brain

Zhenghong Li, Jiaxiang Ren, Zhilin Zou, Kalyan Garigapati, Congwu Du, Yingtian Pan, and Haibin Ling

Stony Brook University, USA  
 {zhenghli,hling}@cs.stonybrook.edu

**Table 1.** Comparison of the influence of the normal B-scan selection threshold.

|               | Normal B-scans |              |              |                   | Corrupted B-scans |              |              |                   |
|---------------|----------------|--------------|--------------|-------------------|-------------------|--------------|--------------|-------------------|
|               | CNR            | MSR          | Dice         | Dice <sub>s</sub> | CNR               | MSR          | Dice         | Dice <sub>s</sub> |
| 1*std         | <b>14.60</b>   | 14.63        | 87.93        | 69.94             | <b>15.33</b>      | 17.12        | 87.17        | 65.04             |
| 1.5*std       | 14.50          | 15.28        | <b>88.18</b> | <b>70.52</b>      | 14.92             | 17.63        | <b>87.64</b> | <b>65.91</b>      |
| 2*std         | 14.56          | 15.62        | 87.80        | 69.77             | 15.04             | 18.22        | 87.30        | 65.29             |
| w/o selection | 14.24          | <b>15.69</b> | 87.67        | 69.51             | 13.75             | <b>19.09</b> | 84.81        | 60.79             |



**Fig. 1.** Comparison of normal B-scan selection thresholds. Stripes are marked by green arrows, and removed stripes are marked by yellow arrows. With reasonable thresholds for normal B-scan selection for training, SOAD can remove the BMA in several consecutive corrupted B-scans. However, without normal B-scan selection, SOAD may not be effective in removing BMA in such cases.

**Table 2.** Comparison of the influence of  $\alpha$  in the weighted loss.


|                 | Normal B-scans |              |              |                   | Corrupted B-scans |              |              |                   |
|-----------------|----------------|--------------|--------------|-------------------|-------------------|--------------|--------------|-------------------|
|                 | CNR            | MSR          | Dice         | Dice <sub>s</sub> | CNR               | MSR          | Dice         | Dice <sub>s</sub> |
| $\alpha = 10$   | 14.31          | <b>15.71</b> | 87.93        | 69.97             | 14.76             | <b>17.94</b> | 87.20        | 65.23             |
| $\alpha = 100$  | <b>14.50</b>   | 15.28        | <b>88.18</b> | <b>70.52</b>      | <b>14.92</b>      | 17.63        | <b>87.64</b> | <b>65.91</b>      |
| $\alpha = 1000$ | 14.27          | 14.67        | 88.00        | 70.31             | 14.76             | 16.70        | 87.45        | 65.46             |

**Table 3.** Comparison of the influence of  $\gamma$  in the weighted loss.

|                | Normal B-scans |              |              |                   | Corrupted B-scans |              |              |                   |
|----------------|----------------|--------------|--------------|-------------------|-------------------|--------------|--------------|-------------------|
|                | CNR            | MSR          | Dice         | Dice <sub>s</sub> | CNR               | MSR          | Dice         | Dice <sub>s</sub> |
| $\gamma = 1$   | 14.06          | 14.32        | 87.98        | 70.07             | 14.36             | 16.64        | 86.90        | 64.43             |
| $\gamma = 1/2$ | 14.14          | 14.57        | 88.09        | 70.47             | 14.62             | 16.70        | 87.41        | 65.54             |
| $\gamma = 1/3$ | <b>14.50</b>   | <b>15.28</b> | <b>88.18</b> | <b>70.52</b>      | <b>14.92</b>      | <b>17.63</b> | 87.64        | 65.91             |
| $\gamma = 1/5$ | 14.27          | 14.68        | 87.85        | 69.82             | 14.75             | 16.69        | <b>87.66</b> | <b>66.10</b>      |

**Table 4.** Comparison of the influence of the window size  $T$ .

|         | Normal B-scans |              |              |                   | Corrupted B-scans |              |              |                   |
|---------|----------------|--------------|--------------|-------------------|-------------------|--------------|--------------|-------------------|
|         | CNR            | MSR          | Dice         | Dice <sub>s</sub> | CNR               | MSR          | Dice         | Dice <sub>s</sub> |
| $T = 5$ | 13.87          | 14.24        | 87.55        | 69.10             | 14.32             | 16.52        | 86.98        | 64.16             |
| $T = 7$ | <b>14.50</b>   | <b>15.28</b> | <b>88.18</b> | <b>70.52</b>      | <b>14.92</b>      | <b>17.63</b> | <b>87.64</b> | <b>65.91</b>      |
| $T = 9$ | 13.91          | 14.62        | 87.83        | 69.86             | 14.36             | 16.79        | 87.33        | 65.86             |

**Fig. 2.** 3D visualization of the raw volume and SOAD (S) results. SOAD can effectively remove the noise in OCTA volumes and improve the 3D visual quality of the volume.